

GM11

Gingiva Mask

Gingiva Mask GM11 is an excellent choice for replicating gingiva contours for digital implant cases. GM11 is a flexible material that can assist in planning the emergence profile for digital implant restorations.

→ Flexible, gingiva color



■ Physical Properties – Liquid Status

Color	Transparent Pink
Viscosity @25°C cps	1179
Wavelength	405nm

■ Mechanical Properties – Post Cured Parts

Property	Value	Test Method
Hardness Shore A	70	ISO 868
Tensile Modulus MPa	N/A	ISO 527
Tensile Strength MPa	N/A	ISO 527
Elongation at break %	N/A	ISO 527
Flexural Modulus MPa	N/A	ISO 178
Flexural Strength MPa	N/A	ISO 178

Post Processing	Value
Cleaning time mins	6
Recommended curing time mins	3 mins per side